MUTSUMI

Hiro®

BFE
≥95%

Disposable Surgical Masks

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(EN) Disposable Surgical Masks

- Helps reduce exposure by body fluids and blood;
- Minimizes medical personnel and patient contamination to exhaled microorganisms;
- 3 Layered design, 3 ply: hydrophobic non-woven outer layer, meltblown filter middle layer and absorbent spunbond inner layer;
- With elastic earloop;
- Adjustable nose bridge for enhanced comfort and fit;
- Bacterial filtration efficiency BFE_≥95%;
- Differential Pressure <40Pa/cm²;
- EN ISO 14683:2019 + AC:2019 -Type I.

(EN) How to use:

- Clean your hands with soap and water or hand sanitizer before touching the mask.
- Remove a mask from the box by the earloops.
- Determine which side of the mask is the top. The side of the mask that has a stiff bendable edge is the top and is meant to mold to the shape of your nose.
- The colored side of the mask is the front and should face away from you, while the white side touches your face.
- Hold the mask by the ear loops. Place a loop around each ear.
- Mold or pinch the stiff edge to the shape of your nose.
- Pull the bottom of the mask over your mouth and chin.

(EN) How to remove a face mask:

- Clean your hands with soap and water or hand sanitizer before touching the mask. Avoid touching the front of the mask. The front of the mask is contaminated.
- Only touch the ear loops.
- Hold both of the ear loops and gently lift and remove the mask.
- Throw the mask in the trash. Clean your hands with soap and water or hand sanitizer.

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(EN) Store in a clean, cool, dry place away from direct sunlight.

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